



Sanquin

# TMB substrate for PeliSPOT

## Specification sheet

<b>Art.no</b>	M2521
<b>Application</b>	ELISPOT
<b>Packing</b>	Each bottle contains 18 ml 1.4 mM 3,3',5,5'-Tetramethylbenzidin / 2.4 mM hydrogen peroxide, pH 5.0
<b>Preservative</b>	none.
<b>Storage</b>	TMB-substrate should be stored at 2-8°C in the dark. Do not use this product beyond the stated expiration date.
<b>Directions for use</b>	Ready for use.

Reagents

Plesmanlaan 125

1066 CX Amsterdam

The Netherlands

Phone: +31 20 512 3599

Fax: +31 20 512 3570

E-mail: [reagents@sanquin.nl](mailto:reagents@sanquin.nl)

Website: [www.sanquinreagents.com](http://www.sanquinreagents.com)

*For research only! Not for human use!*

M2521\01 1808081413